



SageHLS Cassette Kit Workflow Guide

High Molecular Weight DNA Extraction PN# HEX-0004 or HEX-0012

Reagents Supplied by Sage Science			Storage Temp.
4 / 12 ea.	Agarose gel cassettes		RT
20 / 60 ea.	Adhesive Tape Strips		N/A
1 ea.	HLS Lysis Reagent 3% SDS, 10 / 30 ml	A	RT
1 ea.	Enzyme Buffer, 15 / 40 ml	C	4°C
1 ea.	Running Buffer, 40 / 115 ml	E	RT
1 ea.	HLS Lysis Reagent 3% Sarkosyl, 10 / 30 ml	G	RT
1 ea.	HLS Lysis Reagent 1% SDS, 10 / 30 ml	H	RT

Materials Supplied or Prepared by User	Supplier	Cat#
NEBNext® dsDNA Fragmentase®	New England Biolabs	M0348S


Important!

- Prior to day of extraction, schedule availability of cells. Check that cell preparation reagents are ready.
- Maximum recommended cell load per lane for HLS cassettes contains **2.5-10 ug of genomic DNA**, (approximately 1.5e06 human diploid cells).
- Irrespective of input cell load, final input sample volume should be fixed at **70 ul per lane**.

HEX kit Workflow Guide Revision Change Log

Last Rev	New Rev	Date	Page#	Notes
C1	D	12/12/17	12-2	Added Revision Change Log
D	E	6/18/18	12-1	Updated NEB fragmentase part number
D	E	6/18/18	12-7	Revise Workflow File Names
E	F	8/17/18	12-6	Replaced buffer chamber image to show even fill line.
E	F	8/17/18	12-12	Corrected dilution and math for E.coli fragmentase dilution.



HMW DNA Extraction Workflow Summary



A. Prepare cells using Cell Suspension Guide (Appendix A)


B. Prepare the Gel Cassette(s)

1. Remove the gel cassette from the foil bag.
2. **Before removing tape!** – Hold the cassette with top surface almost vertical and the elution modules above the gel channel. Tap the cassette to dislodge any bubbles that are trapped behind the elution wells, or in the elution channels. Repeat if necessary. Allow the bubbles to collect in the electrode channel directly above.




Bubbles in the elution paths can interfere with collection

- 3 Slowly rotate the cassette to allow the bubbles to collect in the upper buffer area. Gently tap if necessary.




Move any bubbles to the upper buffer area.




Aggregate and collect air bubbles in the upper buffer chambers

4. With cassette held at a slight angle to keep bubbles located in the upper buffer chamber, gently place the cassette onto the nest. The upper buffer chamber should be on the left side of the nest.
5. While holding the cassette firmly in place on the nest, grab the tape tab, and pull the tape off at an angle, slowly and firmly. Alternate the pulling angle if the tape resists peeling.




Peel back tape at an angle. Alternate angles if the tape resists peeling.

6. Remove all buffer from all elution wells (set a pipette to 100 μ l to completely empty the wells). Keep pipette tips vertical in the well to avoid damage to the membranes.
7. Taking care not to introduce additional bubbles into the elution modules, add 80 μ l of buffer to all elution wells.



8. Add **Running Buffer** to the upper buffer chamber on each side of the gel cassette until the level is flush with the top of the cassette. An adequate buffer level is important for achieving best results from the SageHLS.



Important! Fill until the buffer level visually reaches the bottom side of the cassette cover.


C. Load the Workflow File

1. Use the following Table as a guide to select the most appropriate Workflow File:


Core Workflow* File Name	Description	Run Time
HMW DNA Extraction	Maximum separation, ultra-HMW fragments in well 2	5:00
HMW High-Pass 50kb elu45	Compression band, 50kb and above, in wells 2 and 3, lower yield	3:30
HMW High-Pass 50kb	Compression band, 50kb and above, in wells 2 and 3	4:20
HMW High-Pass 300kb	3 hour extraction, fragments >300 in well 2	13:00
Non-Core Workflow File Name	Description	Run Time
HMW High-Pass 250kb	3 hour extraction, fragments >250 in well 2	13:00
HMW High-Pass 350kb	3 hour extraction, fragments >350 in well 2	13:00
HMW High-Pass 500kb	3 hour extraction, fragments >500 in well 2	13:00

*Core work flows have been extensively tested with biological samples. Non-core workflows have only been tested with model DNA samples.


2. Go to the Main screen of the SageHLS software:



3. Select the Workflow File from the drop down menu.



- Choose the lanes to be used by clicking the boxes next to the lane numbers and enter sample IDs into the adjacent fields (sample IDs are optional or can be entered later).



Check marks indicate which lanes are active

D. Run the Check Current Test

- Press the "Check Current" button.



- A pop-up window will appear. Press "Start" to begin the Check Current routine.



Press "Start"

- The routine will first test the separation electrodes, then test the elution electrodes, and complete within a few minutes. After a successful test, all boxes will fill with green check marks. Press "Return" to continue




Press "Return"

E. Stage 1: Extraction

1. Use the following Table as a guide to select the most appropriate Lysis Reagent. The extraction step will take 3 hours:


Lysis Reagent Name	Description
HLS Lysis Reagent 3% SDS (A)	Maximum DNA recovery, residual SDS (0.01-0.03%) in eluant
HLS Lysis Reagent 1% SDS (H)	Less DNA recovery (25-30%), less residual SDS (0.001 – 0.009%)
HLS Lysis Reagent 3% Sarkosyl (G)	Less DNA recovery (25-30%), for samples with potassium (alternative to SDS, which co-precipitates in the presence of potassium)

2. Empty all sample and reagent wells. Use caution not to pierce agarose at the bottom of the wells.




3. Load samples in all lanes. Always use 70ul sample loading volume. (Sample wells will not be completely full.)
4. Fill Reagent Wells with HLS Lysis Buffer. Fill, but do not overfill! Leave a concave meniscus to prevent contact with sealing tape in next step. Approximate volume needed will be 220-230 ul.

- Seal reagent, sample, and elution ports with supplied tape without occluding the electrode ports. Press tape firmly around edges of the ports.



- Close the lid and press "Run Workflow". The Extraction step will take 3 hours of unattended operation.



F. Stage 2: Treatment Stage

- At the end of the Extraction Stage/Step, the SageHLS will pause on the first step of the Treatment Stage and a pop-up window will appear with user instructions.

The screenshot displays the SageHLS software interface. At the top, there are tabs for 'Main Screen', 'Workflow Editor', 'Log Review', and 'System Options'. The 'Workflow Editor' tab is active, showing a 'Workflow File' field with 'HLS_HMW DNA Extraction' and a 'Next Configuration' section with several input fields and dropdown menus. Below this, there are status indicators for 'Voltage, V' (0.0), 'LF Status', 'Lock Status', and 'System Status'. A 'Step Table' is visible on the left, listing steps from 1-001 to 3-004. Step 2-001 is highlighted in yellow. A graph on the right shows 'Electrode Current and Sample Well Temperature' over time. A blue pop-up window is overlaid on the graph, containing the following text: 'Remove tape seal. Replace sample well contents with 80ul of reaction mix. Replace reagent well contents with 230ul Enzyme Buffer (B). When done, close the lid and click OK to continue.' Below the text is an 'OK [Enter]' button. Red arrows point from the text labels below to the corresponding elements in the screenshot.

Instrument pauses at 2-001 (stage 2, step 1)

User instructions

Remove tape seal. Replace sample well contents with 80ul of reaction mix. Replace reagent well contents with 230ul Enzyme Buffer (B)
When done, close the lid and click OK to continue.

OK [Enter]



Important! The instrument will remain paused while the user action pop-up window is displayed. It is resumed when the “OK” button is pressed. If the instrument is inadvertently resumed, press the “Pause Workflow” button in the Command Menu to re-pause the instrument and continue with the manual user action.

2. Prepare the Enzyme Reaction Mix:

- a. Remove NEB Fragmentase from the freezer, briefly vortex (1s) to mix
- b. Dilute the NEB Fragmentase (NF) with Enzyme Buffer (C) as follows:


Cell Type	Dil.Factor	Fragmentase	Procedure
E.coli	1:600	0.01 μ l / reaction	<ol style="list-style-type: none"> i. add 2μl of NF to 28μl of Enzyme Buffer C, vortex to mix ii. add 20μl of dilution to 780μl of Enzyme Buffer C, vortex to mix
White Blood Cells	1:400	0.02 μ l / reaction	i. add 2 μ l of NF to 798 μ l of EB, vortex to mix



Important! Fragmentase Enzyme Mix should be used within minutes of preparation. It can be prepared at the end of the extraction step and kept on ice. Preparing the mix within 15 minutes of use is recommended.


3. Open the lid and carefully remove the adhesive tape(s). To remove the tapes, grab the tab in right upper corner and peel diagonally with a slow smooth motion.

grab the tab in right upper corner and peel diagonally with a slow smooth motion



Important! Pulling the tape in a diagonal fashion prevents liquid transfer between adjacent elution ports and transfer between the sample/reagent ports and the elution ports.

- Remove all contents from the Reagent wells and Sample wells in the cassettes to be run. The well volumes are 270 μ l and 85 μ l, respectively.
- Add 70 μ l of the Fragmentase Enzyme Mix to the Sample well.
- Add 230 μ l Enzyme Buffer (C) to the Reagent well.
- Close the lid (**do not re-seal the wells with tape**).
- Press “OK” in the pop-up window to resume the workflow.



Press “OK” to resume

- The enzymatic treatment will take **30 minutes**.

- 10, At the end of **30 minutes** the SageHLS will pause, and a pop-up window with user instructions will appear.

Instrument pauses at 2-003 (stage 2, step 3)

User instructions

Stop Reaction: Add 0.3M EDTA to sample well.
Replace reagent well contents with 230ul of Lysis Buffer (A).
When done, close the lid and click OK to continue.

OK [Enter]




Important! The instrument will remain paused while the user action pop-up window is displayed. It is resumed when the “OK” button is pressed. If the instrument is inadvertently resumed, press the “Pause Workflow” button in the Command Menu to re-pause the instrument and continue with the manual user action.

11. Open the lid and remove the contents of the Reagent well.
12. Replace the Reagent well contents with Lysis Reagent (**A**, **G**, or **H**), ~230 ul.
13. Close the lid and **re-seal the cassette wells with tape.**

Tabs for removing tape seals should be on same side as elution modules.


Be sure that all Reagent, Sample, and Elution ports are completely covered.



Seal tape by pressing firmly on the edges of the ports.

Lower edge of tape must not occlude lower electrode port


15. Close the lid and press "OK" in the pop-up window to resume the workflow.



Press "OK" to resume

G. Stage 3: Collection Stage


1. The Collection Stage will require several hours of unattended operations. Users should note the time remaining, after which the samples can be collected.




Time remaining to the end of the run

2. After the run is complete, open the lid and remove the sealing tape from the cassette(s).

grab the tab in right upper corner and peel diagonally with a slow smooth motion



- SageHLS workflows are designed to collect the target range (indicated by target size in HLS-CATCH or highest molecular weight DNA in HMW workflows) in well number 2. This does not guarantee that all targets will be in that well, and users may have multiple size targets. It is highly recommended that the adjacent wells, 1 and 3, are sampled for target DNA at the least.



Workflows are designed to collect most targets in well #2

- Using a wide-bore pipette tip, remove the contents of the elution modules.



Important! Pipette as slowly as possible to avoid shearing the HMW DNA. Use of an electronic pipettor at low speed settings may be helpful. There should be 70-80 ul of liquid in each module.

- Extremely HMW DNA will be very inhomogeneously distributed in the elution product. To quantify, we recommend Qubit assays using at least three 1 ul aliquots from different locations within the tube. Average the three readings. A high average value with a high CV is diagnostic of very HMW DNA.



For Qubit assays, using at least three 1 ul aliquots from different locations within the tube