

BAKER

Environments For Science™

ReCO₂ver™

Rapid Recovery Incubator

ABSOLUTE PRECISION.
MAXIMUM PROTECTION.

PRECISION



CONTROL OVER TEMPERATURE,
CARBON DIOXIDE & RELATIVE HUMIDITY

FAST



RECOVERY OF CELL CULTURE
CONDITIONS

MAXIMUM



PROTECTION FOR CELLS

BETTER THAN



ISO CLASS 4 (CLASS 10) AIR
CLEANLINESS

In all scientific research applications, each cell type will benefit from keeping proper physiological cell growth conditions.

Cell culture incubators attempt to mimic physiological conditions, but atmospheric conditions within the laboratory and the environmental parameters incubators attempt to control vary widely, and are lost with each door opening taking a long time get back to your ideal parameters. This will impact the integrity of your work.



There are several factors an incubator should control in order to help you achieve optimal cell growth conditions. Precision in temperature, gas and humidity as well as the rapid recovery of those conditions after door openings are crucial for ensuring cells are exposed to a constant environment required for their well-being. These parameters can change very quickly once a door is opened, and it can take a long time for *in vivo*-like conditions to be achieved again.

Key Benefits



- The vertical, uni-directional downward airflow with full-face HEPA filter delivers better-than-ISO Class 4 (Class 10) clean air to the chamber to provide unparalleled protection of your work.
- Proprietary IntelliCELL™ P.I.D. control algorithm offers true active humidity control - giving users complete freedom to precisely define three variables critical to cell growth: temperature, CO₂ and relative humidity.
- Superior stability, uniformity and recovery times - mean that even with a full load, cell cultures grow consistently from shelf to shelf, for dependably high-quality, robust cells on every plate.
- Outstanding condensation control virtually eliminates wet spots where contaminants can grow and spread, even at humidity levels above 90%.
- Fogless interior door with heated frame provides a crystal clear view to every shelf, reducing the need for door openings.
- Ultrasonic humidity delivery system lets you eliminate the use of a water pan, and the risk of contaminants that typically go with it.
- Intuitive color touchscreen controls give users full control over incubator functionality.
- H₂O₂ total biodecontamination (optional or on board for ReCO₂ver™ Plus) in addition to UV light, kills even resistant contaminating microbes in approximately four hours (vs. up to 12 hours for other decontamination methods).

Accessories

- Rolling Cart
- Stacking Kit
- CondoCell® Starter & Add-on Kits
- Additional Shelves
(7 additional, 12 max)
- Two-Stage CO₂ Gas Tank Regulator
(adjustable to max of 15 psi)
- Portable CO₂ Meter
& Portable pH Meter
- CO₂ Analyzer Fyrite® Kit
- Refill CO₂ Analyzer Fluid Kit
- Seismic Restraints
- Preventive Maintenance Kit
- Dry In-Line Filter



Specification


		ReCO ₂ ver™		ReCO ₂ ver™ Plus	
Performance Specifications					
Model Number		REC 602		REC 602 -Plus	
Humidity Delivery		Ultrasonic (Nebulizer)			
Heat Type		Direct Heat			
External Dimensions		26.9" W x 26.5" D x 40.5" H (683mm W x 673mm D x 1029mm H)			
Internal Dimensions		19.9" W x 21.9" D x 33.5" H (505mm W x 556mm D x 851mm H)			
Usable Interior Volume / Capacity		5.92 ft ³ (0.168 m ³) / 275 90mm plates			
Usable Shelf Area / Max Shelves		18.6" W x 20.1" D (472mm W x 511mm D) / 12 shelves (4 standard)			
Weight		250 lbs. (114 kg.) / Shipping Weight 310 lbs. (141 kg)			
Environmental Performance					
Temperature Control Range		13°F (7.0°C) above ambient - 131°F (55.0°C) (increments of 0.2°F / 0.1°C)			
Temperature Accuracy		+/-0.1°C			
Temperature Uniformity		+/- 0.25°C			
Temperature Recovery Time		Within 1°C of set point (37°C) in 6 minutes			
CO ₂ Control Range		0.0% - 20.0% (increments of 0.1%)			
CO ₂ Accuracy		+/-0.15%			
CO ₂ Recovery Time		Within 0.2% of set point (5%) in 5 minutes			
Contamination Control					
Relative Humidity Control Range		Two operating modes: Up to 90%; 90%-95%			
Relative Humidity Accuracy		+/-3%			
Relative Humidity Recovery Time		Within 3% of set point (90%) in 4 minutes			
Contamination Prevention		Expansive HEPA filter providing better-than-ISO Class 4 (Class 10) air cleanliness			
Biodecontamination Methods		UV Light		UV Light & H ₂ O ₂ (Standard)	
		ReCO ₂ ver™		ReCO ₂ ver™ Plus	
Features		Standard	Optional	Standard	Optional
HEPA Filtration (better than ISO Class 4)		✓		✓	
Active Humidity Control		✓		✓	
Active Temperature Control		✓		✓	
Active CO ₂ Control		✓		✓	
Alarm Contacts		✓		✓	
UV Light		✓		✓	
Automatic Water Refill System		✓		✓	
Direct Heat – Rapid recovery of temperature back to set point		✓		✓	
H ₂ O ₂ Biodecontamination			✓	✓	
Copper Interior Components			✓		✓
CO ₂ Supply Automatic Switchover			✓		✓
RS232/485 or 4-20mA			✓		✓

Additional equipment compatibility




CondoCell®

Portable device that allows transport from/to the incubator while maintaining the cell culture environment.



Etaluma Lumascope

Live cell imaging.



PhO₂x Box

Cell Culture Chamber with CO₂ and O₂ control.



Resources



Read how the air filtration system inside ReCO₂ver™ performs against a leading competitor.



Uniformity and Consistency – Eliminating the Edge Effect. Download the white paper.



Active humidity and condensation control, even at >90% RH. Download the white paper.



Learn about our H₂O₂ biodecontamination protocol.

BAKER

Environments For Science™

For North American Inquiries:

bakerco@bakerco.com

+1 (800) 992-2537

www.bakerco.com

UK / EMEA and Global Inquiries:

sales@cleanair.eu.com

+31 (0)85 90 22 500

www.cleanair.eu.com