

Technical Support Bulletin

UranylLess classic and negative staining protocols:

[#29-005030](#), [#29-005040](#)


UranylLess – Protocol for classic contrast

This protocol is used for double staining with UranylLess and lead citrate on ultrathin sections.

The protocol is adapted for biological samples that have been fixed with glutaraldehyde, osmium or ruthenium and embedded in an epoxy type resin (i.e. Epon, Araldite, Spurr) or an acrylic type resin (LRWhite, HM20).

Staining protocol:


1. Place a droplet of UranylLess on parafilm or a hydrophobic slide
2. Place the TEM grid on the UranylLess droplet for 1-2 minutes
3. Blot the grid with filter paper
4. Then wash the grid in distilled water
5. Let the grid with sample dry



**Place a droplet of
UranylLess on parafilm**



**Stain section
1-2 minutes**



Blot with filter paper



**Rinse in distilled water
and let sample dry**

After drying, use the lead citrate protocol according to Reynolds (1963).

1. Place the grid on a 3% lead citrate droplet for 1 minute
2. Blot the grid with filter paper
3. Rinse the grid with distilled water
4. Let the grid dry

Technical tips:

- UranylLess is not air or light sensitive (unlike uranyl acetate)
- After lead citrate, rinse immediately in freshly prepared distilled water or wash with a solution containing 0.01N of NaOH
- If there is any precipitate in the solution, filter it prior to use

Technical Support Bulletin


- If solution was refrigerated, allow solution to return to room temperature before use
- Do not keep lead citrate or UranyLess refrigerated

UranyLess – Protocol for negative staining

This protocol is adapted for negative staining with UranyLess. It allows characterisation of isolated particles or morphology of bacteria, viruses, protein, liposomes, exomes, nanoparticles etc.

Staining protocol:

1. Place a droplet of the sample solution (~10ul) and a droplet of UranyLess on a piece of parafilm or any other hydrophobic carrier
2. Place a formvar coated grid on the sample solution for about 1 minute; use fine tweezers.
3. Blot the grid with filter paper
4. Place the grid on the UranyLess droplet for about 1 minute
5. Blot the grid with filter paper
6. Let dry for 5 minutes and observe with a TEM



Place a sample droplet

Place grid on sample droplet

Blot with filter paper

Place the grid on UranyLess

Blot with filterpaper

Technical Tip:

- If the staining is too intense, wash with water for 1 minute.

Technical Support Bulletin

TSB 28-005030 UranyLess classic and negative staining protocols 10-15-2016 Rev.1